



8.

Spectrum Spectrogram Map

Spectrum Info Spectrum ROIs Spectrogram ROIs Real Time Detectors Isotopics Console System Log

Clear ROIs View Full Info CPS Export to CSV Save Reference ROIs ☐ Background Subtracted

	ROI	Spectrum	Lower	Upper	Peak Center	FWHM	Peak Area	ROI Counts	N Counts	B Counts	Nuclide	Photo Peak	Peak D
1	Ba-133, 303 keV	Ba_Cs_Co	193 keV	297 keV	Fit Failed	None	None	2 341 784	None	None	None	None	None
2	Cs-137, 662 keV	Ba_Cs_Co	584 keV	738 keV	661.07 ± 0.15 keV	58.09 ± 0.44 keV	164 940 ± 1 619	481 990	164 502 ± 1 912	317 487 ± 1 782	Cs-137	661.6553 keV	-0.589 keV
3	Co-60, 1173 keV	Ba_Cs_Co	1093 keV	1249 keV	1174.18 ± 1.97 keV	70.02 ± 6.44 keV	46 494 ± 8 685	121 467	46 740 ± 9 428	74 726 ± 9 421	Co-60	1173.228 keV	0.954 keV
4	ROI_8	Ba_Cs_Co	1244 keV	1451 keV	1332.64 ± 2.05 keV	88.31 ± 10.69 keV	58 495 ± 14 178	76 335	58 449 ± 14 053	17 885 ± 14 050	Co-60	1332.492 keV	0.148 keV