

Reset Zoom

Unlock y-axis

Log


CPS

Cursor

Add ROI

No Background

Cursor



3.

ROI Editor

ROI Name	ROI_2
Lower Bound	337,00 keV
Higher Bound	457,00 keV
Peak Function	Gaussian
Background	Quadratic
Edge Channels	3
Nuclide	None
Photopeak	None

Auto Match Nuclide

☒ Allow merging

☒ Movable

☐ Use Poisson Weights

Delete Cancel OK

Spectrum

Spectrogram

Map